


Polyomaviruses BK and JC PCR reagents

Overview:.....	1
Products.....	2
BK Master mix FAM-BHQ1 MM5700 0.175ml \$32.....	2
BK FAM-BHQ1 Primer-probe PP5700 0.055ml \$80.....	2
JC Master mix FAM-BHQ1 MM5800 0.175ml \$32	2
JC FAM-BHQ1 Primer-probe PP5800 0.055ml \$80.....	2
BK Plasmid PLAS5700 0.25ml \$35	2
JC Plasmid PLAS5800 0.25ml \$35	2
Detection of Polyomavirus DNA:.....	3
Thermal cycle conditions for PCR reactions on RotorGene*.....	3
Thermal cycle conditions for PCR reactions on LightCycler.....	3
Use of the Polyomavirus plasmids:.....	3
Attostar Master mix MM5700 or MM5800 in RotorGene*- PCR 24 µl reactions	4
Attostar Master mix MM5700 or MM5800 in LightCycler- PCR 24 µl reactions.....	4
Attostar Primers/probe PP5700 or PP5800 in RotorGene-PCR 20µl reactions	5
Attostar Primers/probe PP5700 or PP5800 in LightCycler-PCR 20 µl reactions.....	5
Polyomavirus and T4 multiplex reaction (PP5700 or PP5800 and PP160).....	6
<i>RotorGene Multiplex PCR 20µl reactions</i>	6
Brief procedure for use of T4 as extraction and amplification control:.....	6
Typical standard curve figures:.....	6
Reference:	7

Overview:



Polyomaviruses are ubiquitous dsDNA human and animal viruses. Most adults have been infected and harbor latent virus with no serious consequences. However immune suppression may result in serious diseases including JC associated progressive multifocal leukoencephalopathy (PML) and BK associated graft failure in renal transplant patients and systemic, lung, liver, and brain disease in other immune-compromised persons. Wu and Ki are newly described human respiratory polyoma viruses.

Polyomaviruses may infect many species including: monkeys, mouse, cows, rabbits, birds, and others. Infections of species other than the natural host may cause tumors, i.e. SV40, a monkey virus, may cause sarcomas and leukemia in hamsters.

Dr Sara Stewart 1906-1976 (Georgetown University's first medical doctor) and Dr Bernice Eddy 1903-1989 were the first to describe polyomaviruses (1). They showed the viruses to cause tumors in many animals.

Master mix, primers, probes, and buffer are available for amplification and detection of the BK and JC. The BK primers are specific for a 140 bp portion of coat protein VP2. This amplification reagent detects BK not JC. The BK detection probe is a molecular beacon. The JC primers are specific for a 140 bp portion of coat protein VP2. The JC probe detects JC not BK. The BK and JC detection probes are molecular beacons. (2)

Products

Products	Catalog No.	Quantity	Price
BK Master mix FAM-BHQ1 <i>Store at -20C. 8 reactions</i> <i>Typical use: Add Taq and DNA, then thermal cycle.</i> <i>Master mix contains buffer, dNTPs, primers, and probe.</i> <i>Detection at 510nm.</i>	MM5700	0.175ml	\$32
BK FAM-BHQ1 Primer-probe <i>Store at -20C. 20 reactions</i> <i>Typical use: Add 2X master mix and DNA, then thermal cycle.</i> <i>Attostar reagent contains primers and probe.</i> <i>Detection at 510nm.</i>	PP5700	0.055ml	\$80
JC Master mix FAM-BHQ1 <i>Store at -20C. 8 reactions</i> <i>Typical use: Add Taq and DNA, then thermal cycle.</i> <i>Master mix contains buffer, dNTPs, primers, and probe.</i> <i>Detection at 510nm.</i>	MM5800	0.175ml	\$32
JC FAM-BHQ1 Primer-probe <i>Store at -20C. 20 reactions</i> <i>Typical use: Add 2X master mix and DNA, then thermal cycle.</i> <i>Attostar reagent contains primers and probe.</i> <i>Detection at 510nm.</i>	PP5800	0.055ml	\$80
BK Plasmid <i>Store at -20C.</i> <i>Typical use: make serial 10 fold dilutions in TE for standard curve, diluting 5 ul into 45 ul water.</i>	PLAS5700	0.25ml	\$35
JC Plasmid <i>Store at -20C.</i> <i>Typical use: make serial 10 fold dilutions in TE for standard curve, diluting 5 ul into 45 ul water.</i>	PLAS5800	0.25ml	\$35

Detection of Polyomavirus DNA:

.....

Thermal cycle conditions for PCR reactions on RotorGene*

*Similar cycle conditions and reaction volumes may be used on many other thermal cyclers.

- 95C 600 sec (activation for AmpliTaq Gold polymerase)
- 40 cycles
 - 95C 15 seconds
 - 60C 30 seconds RotorGene Channel Setup FAM/Sybr, Cy5; Gain 7
 - 72C 30 seconds

FAM/Sybr has a source of 470nm and Detector 510nm (LightCycler use F1)
Cy5 has a source of 625nm and Detector 660hp nm
Quasar 670 has the same fluorescent absorption and emission as Cy5.

.....

Thermal cycle conditions for PCR reactions on LightCycler

- 95C 600 sec (activation for AmpliTaq Gold polymerase)
 - 40 cycles
 - 95C 15 seconds
 - 60C 30 seconds acquire fluorescent signal on F1 gain =1
 - 72C 30 seconds
 - 40C 30 seconds cool
-

Use of the Polyomavirus plasmids:

Dilute the plasmid in TE to prepare a standard curve. Common dilutions would be 10-fold from 200 to 0.02pg/ml. The 0.02 pg/ml plasmid dilution contains 20 copies of plasmid in 2 µl

.....

Attostar Master mix MM5700 or MM5800 in RotorGene*- PCR 24 μ l reactions

*Similar cycle conditions and reaction volumes may be used on many other thermal cyclers.

Reagents needed for 24 μ l final volume reactions:								
Reaction tube number	1	2	3	4	5	6	7	8
Master mix MM5700 or MM5800	19.6	39.2	58.8	78.4	98	117.6	137.2	156.8 μ l
AmpliTaq Gold polymerase	0.4	0.8	1.2	1.6	2	2.4	2.8	3.2 μ l
Dispense 20 μ l /tube								
Add 4 μ l DNA/tube								

Attostar Master mix MM5700 or MM5800 in LightCycler- PCR 24 μ l reactions

Reagents needed for 24 μ l final volume reactions:								
	1	2	3	4	5	6	7	8
Master mix MM5700 or MM5800	19.6	39.2	58.8	78.4	98	117.6	137.2	156.8 μ l
AmpliTaq Gold polymerase	0.4	0.8	1.2	1.6	2	2.4	2.8	3.2 μ l
BSA 1 mg/ml	1	2	3	4	5	6	7	8 μ l
Dispense 21 μ l / capillary								
Add 3 μ l DNA / capillary								

·
·

.....

Attostar Primers/probe PP5700 or PP5800 in RotorGene-PCR 20µl reactions

Reagents needed for 20 µl PCR final reaction tube volumes										
Reaction tube number	1	2	3	4	5	6	7	8	9	10
Attostar Primer-Probe PP5700 or PP5800	2	4	6	8	10	12	14	16	18	20 µl
2 x master mix	12	24	36	48	60	72	84	96	108	120 µl
Dispense 14 µl / reaction tube										
Add 6 µl DNA / reaction tube										

Reaction tube number	11	12	13	14	15	16	17	18	19	20
Attostar Primer-Probe PP5700 or PP5800	22	24	26	28	30	32	34	36	38	40 µl
2 x master mix	132	144	156	168	180	192	204	216	228	240 µl
Dispense 14 µl / reaction tube										
Add 6 µl DNA / reaction tube										

Attostar Primers/probe PP5700 or PP5800 in LightCycler-PCR 20 µl reactions

LightCycler Reagents needed for 20 µl PCR final reaction tube volumes										
Reaction tube number	1	2	3	4	5	6	7	8	9	10
Attostar Primer-Probe PP5700 or PP5800	2	4	6	8	10	12	14	16	18	20 µl
2 x master mix	12	24	36	48	60	72	84	96	108	120 µl
BSA 1 mg/ml	1	2	3	4	5	6	7	8	9	10 µl
Dispense 15 µl / reaction tube										
Add 5 µl DNA / reaction tube										

Reaction tube number	11	12	13	14	15	16	17	18	19	20
Attostar Primer-Probe PP5700 or PP5800	22	24	26	28	30	32	34	36	38	40 µl
2 x master mix	132	144	156	168	180	192	204	216	228	240 µl
BSA 1 mg/ml	11	12	13	14	15	16	17	18	19	20 µl
Dispense 15 µl / reaction tube										
Add 5 µl DNA / reaction tube										

Polyomavirus and T4 multiplex reaction (PP5700 or PP5800 and PP160)

RotorGene Multiplex PCR 20µl reactions



Brief procedure for use of T4 as extraction and amplification control:

- Add 5µl bacteriophage to the sample. Proceed with DNA extraction. Dilutions of the bacteriophage may be made to give a final PCR Ct value that is about 30-35. At this dilution, the phage is more sensitive, i.e. more likely, to detect a poor extraction or the presence of PCR inhibitors in the reaction. Please refer to the PP160 product literature.

	Multiplex Reagents needed for 20 µl PCR final reaction tube volumes										
Reaction tube number	1	2	3	4	5	6	7	8	9	10	
T4 Primer Probe PP160	2	4	6	8	10	12	14	16	18	20	20µl
Polyoma Primer-Probe PP5700 or PP5800	2	4	6	8	10	12	14	16	18	20	20µl
2 x master mix	10	20	30	40	50	60	70	80	90	100	100µl
Dispense 14 µl / reaction tube											
Add 6 µl DNA / reaction tube											

	11	12	13	14	15	16	17	18	19	20	
T4 Primer Probe PP160	22	24	26	28	30	32	34	36	38	40	40µl
Polyomavirus Primer-Probe PP5700 or PP5800	22	24	26	28	30	32	34	36	38	40	40µl
2 x master mix	110	120	130	140	150	160	170	180	190	200	200µl
Dispense 14 µl / reaction tube											
Add 6 µl DNA / reaction tube											

Typical standard curve figures:



.....

Reference:

- 1) Eddy, BE and Stewart SE, Characteristics of the SE Polyoma Virus, Am J Public Health Nations Health 1959 November, 49: 1491-1492.
- 2) <http://www.molecular-beacons.org/Introduction.html>

For research use only. Not intended for any animal or human therapeutic or diagnostic use. This product is sold under license from the Public Health Research Institute. It may be used under PHRI Patent Rights only for the purchaser's research and development activities. 'Black Hole Quencher,' 'CAL Fluor,' 'Pulsar' and 'Quasar' are trademarks of and licensed by Biosearch Technologies, Inc., Novato, CA. The BHQ, CAL Fluor, Pulsar and Quasar dye technology is the subject of existing or pending patents including US Patent No. 7,019,129 and is licensed and sold under agreement with Biosearch Technologies, Inc.

For technical support, contact Attostar@Attostar.com